



LIFECARE SOFTWARE SOLUTIONS, INC.
Headquarters: P.O. Box 173 • Budd Lake, NJ • 07828
Phone: 908-852-3258 • Fax: 908-852-2255

IL Branch: 1365 Wiley Road, Unit 154 • Schaumburg, IL • 60173
Phone: 847-885-8155 • Fax: 847-885-8156
Corporate Email: info@lifecare.net • On the Web: www.lifecare.net

LifeCare Plus® version 2.2

Installation Instruction for UNIX Client

Before Installing

- Verify that you are running LifeCare Plus version 2.1 and ProvideX version 5.14. These version numbers are located at the bottom of the LifeCare Plus Main Menu. ***If you are running an earlier version of LifeCare Plus or ProvideX, Please call LifeCare Software Solutions Support at (847) 885-8155 to discuss your installation options.***
- Verify that you have a **current** backup of your LifeCare directory.
- Be sure you read and understand the Installation Instructions and Release Notes.
- Make sure all users have exited LifeCare before installing.
- Right-click on the LifeCare Plus icon and left-click "Properties". Write down the path in the "Start In" field. For example: **P : \PVX**
- Write down the server IP address in the "Target" field. For example: 192.168.1.144
- Look for "-DR=/u/pvx" or something similar in your "Target" field. Write down what follows "-DR=". For example: /u/pvx
- This is where ProvideX is installed on the server.

Install:

- Double-click the installation program you downloaded to begin the installation procedure.
- Click "Next" when you have finished reading the screen, or click "Cancel" to stop the installation.
- Confirm that the path in the installation folder matches the path in your "Start In" field from above. If not, change it accordingly, or click "Cancel" to stop the installation.
- You will receive a message warning you that the folder already exists. This is expected. Click "Yes" to proceed or "No" to cancel the procedure.
- Click "Next" to start copying files, or click "Cancel" to stop the installation. **Do not click "Cancel" once the copy-files procedure has begun.**
- Click "Ok" to complete the installation.

Grant file permissions to the pvx directory within Unix:

- **At this point you will need the root password for your server - if you don't already have the root password, consult your system administrator before proceeding to the next step.**
- Click the Windows "Start" button. Click "Run".
- Type **telnet <IP ADDRESS>** where <IP ADDRESS> is the server IP address you obtained above. For example: **telnet 192.168.1.144**
- Login to your server as *root*.
- Change to the directory preceding where ProvideX is installed. For example, if ProvideX is installed in /u/pvx you would change to /u using the following command: **cd /u**
- Execute the following command:
- **chmod -R 777 pvx**
- Type **exit** to close your telnet session.



LIFECARE SOFTWARE SOLUTIONS, INC.
Headquarters: P.O. Box 173 • Budd Lake, NJ • 07828
Phone: 908-852-3258 • Fax: 908-852-2255

IL Branch: 1365 Wiley Road, Unit 154 • Schaumburg, IL • 60173
Phone: 847-885-8155 • Fax: 847-885-8156
Corporate Email: info@lifecare.net • On the Web: www.lifecare.net

Update data dictionary files:

- Run LifeCare Plus from the desktop icon.
- Double-click your LifeCare icon to run the software. A series of programs will run to update your current data file structure. Your LifeCare directory will be displayed on the screen called *Location of Database Files*. The default will be printed below the phrase **LifeCare directories**.
There will be at least 4 directories listed, which will all begin with the root LifeCare directory (for example, **P:\PVX\CUSTOM\ P:\PVX\LIB\ P:\PVX\DATA\ P:\PVX\LIBRARY** - where **P:\PVX** is the root LifeCare directory). In standard installs, the information listed will be correct. If the correct LifeCare directory is listed, click **OK**. If not, enter the proper directory in the space provided under **Alternate Location**. You will need to enter the full path of any additional directories (i.e., **P:\PVX\DATA2**), and click **OK**. The database update may take some time, depending on the number of facilities you have installed.
- If you have installed a point release, verify that the version number has updated on the Main Menu. If you have installed a patch, verify that the expected changes are present.